

Item Name	Item No.	Description	Price
pNR-GAL4(DBD)/AR(LBD)	PNR219301	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and AR (AR/NR3C4) LBD domain (ligand-binding domain). [AR, androgen receptor]	Inquiry
pNR-GAL4(DBD)/ERa(LBD)	PNR219302	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and ER alpha (ERA/ESR1/NR3A1) LBD domain (ligand-binding domain). [ESR1, estrogen receptor 1]	Inquiry
pNR-GAL4(DBD)/ERb(LBD)	PNR219303	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and ER beta (ERB/ESR2/NR3A2) LBD domain (ligand-binding domain). [ESR2, estrogen receptor 2]	Inquiry
pNR-GAL4(DBD)/ERRa(LBD)	PNR219304	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and ERR alpha (ERR1/ESRRA/NR3B1) LBD domain (ligand-binding domain). [ESRRA, estrogen related receptor alpha]	Inquiry
pNR-GAL4(DBD)/FXR(LBD)	PNR219305	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and FXR (NR1H4) LBD domain (ligand-binding domain). [NR1H4, nuclear receptor subfamily 1 group H member 4]	Inquiry
pNR-GAL4(DBD)/GR(LBD)	PNR219306	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and GR (NR3C1) LBD domain (ligand-binding domain). [NR3C1, nuclear receptor subfamily 3 group C member 1]	Inquiry
pNR-GAL4(DBD)/LXRa(LBD)	PNR219307	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and LXR alpha (LXRA/N41H3) LBD domain (ligand-binding domain). [NR1H3, nuclear receptor subfamily 1 group H member 3]	Inquiry
pNR-GAL4(DBD)/LXRb(LBD)	PNR219308	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and LXR beta (LXRB/NR1H2) LBD domain (ligand-binding domain). [NR1H2, nuclear receptor subfamily 1 group H member 2]	Inquiry
pNR-GAL4(DBD)/MR(LBD)	PNR219309	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and MR (NR3C2) LBD domain (ligand-binding domain). [NR3C2, nuclear receptor subfamily 3 group C member 2]	Inquiry
pNR-GAL4(DBD)/PPARa(LBD)	PNR219310	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and PPAR alpha (PPARA/NR1C1) LBD domain (ligand-binding domain). [PPARA, peroxisome proliferator activated receptor alpha]	Inquiry

pNR-GAL4(DBD)/PPARd(LBD)	PNR219311	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and PPAR delta (PPARD/NR1C2) LBD domain (ligand-binding domain). [PPARD,peroxisome proliferator activated receptor delta]	Inquiry
pNR-GAL4(DBD)/PPARg(LBD)	PNR219312	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and PPAR gamma (PPARG/NR1C3) LBD domain (ligand-binding domain). [PPARG,peroxisome proliferator activated receptor gamma]	Inquiry
pNR-GAL4(DBD)/PR(LBD)	PNR219313	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and PR (PGR/NR3C3) LBD domain (ligand-binding domain). [PGR,progesterone receptor]	Inquiry
pNR-GAL4(DBD)/RARa(LBD)	PNR219314	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and RAR alpha (RARA/NR1B1) LBD domain (ligand-binding domain). [RARA,retinoic acid receptor alpha]	Inquiry
pNR-GAL4(DBD)/RARb(LBD)	PNR219315	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and RAR beta (RARB/NR1B2) LBD domain (ligand-binding domain). [RARB,retinoic acid receptor beta]	Inquiry
pNR-GAL4(DBD)/RARg(LBD)	PNR219316	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and RAR gamma (RARG/NR1B3) LBD domain (ligand-binding domain). [RARG,retinoic acid receptor gamma]	Inquiry
pNR-GAL4(DBD)/RXRa(LBD)	PNR219317	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and RXR alpha (RXRA/NR2B1) LBD domain (ligand-binding domain). [RXRA,retinoid X receptor alpha]	Inquiry
pNR-GAL4(DBD)/RXRb(LBD)	PNR219318	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and RXR beta (RXRB/NR2B2) LBD domain (ligand-binding domain). [RXRB,retinoid X receptor beta]	Inquiry
pNR-GAL4(DBD)/TRa(LBD)	PNR219319	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and TR alpha (THRA/NR1A1) LBD domain (ligand-binding domain). [THRA,thyroid hormone receptor alpha]	Inquiry
pNR-GAL4(DBD)/TRb(LBD)	PNR219320	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and TR beta (THRB/NR1A2/T3RB) LBD domain (ligand-binding domain). [THRB,thyroid hormone receptor beta]	Inquiry

pNR-GAL4(DBD)/VDR(LBD)	PNR219321	Plasmid expressing the fusion protein of GAL4 DBD (DNA-binding domain) domain and VDR (NR111) LBD domain (ligand-binding domain). [VDR, vitamin D receptor]	Inquiry
------------------------	-----------	---	---------