

Item Name	Item No.	Description	Price
pPRO-CD79B Plasmid	PTS001	Expression plasmid vector containing CD79B promoter sequence for specific reporter gene or GOI expression in B cells	Inquiry
pPRO-CD14 Plasmid	PTS002	Expression plasmid vector containing CD14 promoter sequence for specific reporter gene or GOI expression in Monocytes	Inquiry
pPRO-SPN Plasmid	PTS003	Expression plasmid vector containing SPN promoter sequence for specific reporter gene or GOI expression in Leukocytes, platelets	Inquiry
pPRO-PTPRC Plasmid	PTS004	Expression plasmid vector containing PTPRC promoter sequence for specific reporter gene or GOI expression in Haematopoietic cells	Inquiry
pPRO-CD68 Plasmid	PTS005	Expression plasmid vector containing CD68 promoter sequence for specific reporter gene or GOI expression in Macrophages	Inquiry
pPRO-FLT1 Plasmid	PTS006	Expression plasmid vector containing FLT1 promoter sequence for specific reporter gene or GOI expression in Endothelial cells	Inquiry
pPRO-GFAP Plasmid	PTS007	Expression plasmid vector containing GFAP promoter sequence for specific reporter gene or GOI expression in Astrocytes	Inquiry
pPRO-ICAM2 Plasmid	PTS008	Expression plasmid vector containing ICAM2 promoter sequence for specific reporter gene or GOI expression in Endothelial cells and megakaryocytes	Inquiry
pPRO-SFTPb Plasmid	PTS009	Expression plasmid vector containing SFTPb promoter sequence for specific reporter gene or GOI expression in Lung cells	Inquiry
pPRO-SYN1 Plasmid	PTS010	Expression plasmid vector containing SYN1 promoter sequence for specific reporter gene or GOI expression in Neurons	Inquiry
pPRO-CAMK2G Plasmid	PTS011	Expression plasmid vector containing CAMK2G promoter sequence for specific reporter gene or GOI expression in Brain cells	Inquiry
pPRO-BIRC5 Plasmid	PTS012	Expression plasmid vector containing BIRC5 promoter sequence for specific reporter gene or GOI expression in Cancer cells	Inquiry
pPRO-AFP Plasmid	PTS013	Expression plasmid vector containing AFP promoter sequence for specific reporter gene or GOI expression in Hepatoma cells	Inquiry
pPRO-NPHS1 Plasmid	PTS014	Expression plasmid vector containing NPHS1 promoter sequence for specific reporter gene or GOI expression in Podocytes	Inquiry
pPRO-MB Plasmid	PTS015	Expression plasmid vector containing MB promoter sequence for specific reporter gene or GOI expression in Muscle cells	Inquiry
pPRO-WFDC2 Plasmid	PTS016	Expression plasmid vector containing WFDC2 promoter sequence for specific reporter gene or GOI expression in Ovarian cancer cells	Inquiry
pPRO-CCKAR Plasmid	PTS017	Expression plasmid vector containing CCKAR promoter sequence for specific reporter gene or GOI expression in Pancreatic cancer cells	Inquiry

pPRO-KLK3 Plasmid	PTS018	Expression plasmid vector containing KLK3 promoter sequence for specific reporter gene or GOI expression in Prostate cells	Inquiry
pPRO-GATA2 Plasmid	PTS019	Expression plasmid vector containing GATA2 promoter sequence for specific reporter gene or GOI expression in Stem cells	Inquiry
pPRO-ENG Plasmid	PTS020	Expression plasmid vector containing ENG promoter sequence for specific reporter gene or GOI expression in Endothelial cells	Inquiry
pPRO-FN1 Plasmid	PTS021	Expression plasmid vector containing FN1 promoter sequence for specific reporter gene or GOI expression in Differentiating cells, healing tissues	Inquiry
pPRO-ITGA2B Plasmid	PTS022	Expression plasmid vector containing ITGA2B promoter sequence for specific reporter gene or GOI expression in Megakaryocytes	Inquiry
pPRO-WAS Plasmid	PTS023	Expression plasmid vector containing WAS promoter sequence for specific reporter gene or GOI expression in Hematopoietic cells	Inquiry
pPRO-ERBB2 Plasmid	PTS024	Expression plasmid vector containing ERBB2 promoter sequence for specific reporter gene or GOI expression in Breast & pancreatic cancer	Inquiry
pPRO-PTGS2 Plasmid	PTS025	Expression plasmid vector containing PTGS2 promoter sequence for specific reporter gene or GOI expression in Tumor cells	Inquiry
pPRO-CXCR4 Plasmid	PTS026	Expression plasmid vector containing CXCR4 promoter sequence for specific reporter gene or GOI expression in Tumor cells	Inquiry
pPRO-E2F1 Plasmid	PTS027	Expression plasmid vector containing E2F1 promoter sequence for specific reporter gene or GOI expression in Tumor cells	Inquiry
pPRO-LCP1 Plasmid	PTS028	Expression plasmid vector containing LCP1 promoter sequence for specific reporter gene or GOI expression in Tumor cells	Inquiry
pPRO-MUC1 Plasmid	PTS029	Expression plasmid vector containing MUC1 promoter sequence for specific reporter gene or GOI expression in Carcinoma cells	Inquiry
pPRO-APOA2 Plasmid	PTS030	Expression plasmid vector containing APOA2 promoter sequence for specific reporter gene or GOI expression in Hepatocytes( from hepatocyte progenitors)	Inquiry
pPRO-SERPINA1 Plasmid	PTS031	Expression plasmid vector containing SERPINA1 promoter sequence for specific reporter gene or GOI expression in Hepatocytes( from definitive endoderm stage)	Inquiry
pPRO-CYP3A4 Plasmid	PTS032	Expression plasmid vector containing CYP3A4 promoter sequence for specific reporter gene or GOI expression in Mature hepatocytes	Inquiry
pPRO-INS Plasmid	PTS033	Expression plasmid vector containing INS promoter sequence for specific reporter gene or GOI expression in Pancreatic beta cells (from definitive endoderm stage)	Inquiry
pPRO-IRS2 Plasmid	PTS034	Expression plasmid vector containing IRS2 promoter sequence for specific reporter gene or GOI expression in Pancreatic beta cells	Inquiry
pPRO-PPY Plasmid	PTS035	Expression plasmid vector containing PPY promoter sequence for specific reporter gene or GOI expression in PP pancreatic cell (gamma cells)	Inquiry

pPRO-MYH6 Plasmid	PTS036	Expression plasmid vector containing MYH6 promoter sequence for specific reporter gene or GOI expression in Late differentiation marker of cardiac muscle cells (atrial specificity)	Inquiry
pPRO-MYL2 Plasmid	PTS037	Expression plasmid vector containing MYL2 promoter sequence for specific reporter gene or GOI expression in Late differentiation marker of cardiac muscle cells (ventricular specificity)	Inquiry
pPRO-TNNI3 Plasmid	PTS038	Expression plasmid vector containing TNNI3 promoter sequence for specific reporter gene or GOI expression in Cardiomyocytes (from immature state)	Inquiry
pPRO-NPPA (ANF) Plasmid	PTS039	Expression plasmid vector containing NPPA (ANF) promoter sequence for specific reporter gene or GOI expression in Atrial specificity in adult cells	Inquiry
pPRO-SLC8A1 Plasmid	PTS040	Expression plasmid vector containing SLC8A1 promoter sequence for specific reporter gene or GOI expression in Cardiomyocytes from early development stages	Inquiry
pPRO-SYN1 (hSyn) Plasmid	PTS041	Expression plasmid vector containing SYN1 (hSyn) promoter sequence for specific reporter gene or GOI expression in Neurons	Inquiry
pPRO-GFAP Plasmid	PTS042	Expression plasmid vector containing GFAP promoter sequence for specific reporter gene or GOI expression in Astrocytes	Inquiry
pPRO-INA Plasmid	PTS043	Expression plasmid vector containing INA promoter sequence for specific reporter gene or GOI expression in Neuroprogenitors	Inquiry
pPRO-MOBP Plasmid	PTS044	Expression plasmid vector containing MOBP promoter sequence for specific reporter gene or GOI expression in Oligodendrocytes	Inquiry
pPRO-MBP Plasmid	PTS045	Expression plasmid vector containing MBP promoter sequence for specific reporter gene or GOI expression in Oligodendrocytes	Inquiry
pPRO-TH Plasmid	PTS046	Expression plasmid vector containing TH promoter sequence for specific reporter gene or GOI expression in Dopaminergic neurons	Inquiry
pPRO-FLG Plasmid	PTS047	Expression plasmid vector containing FLG promoter sequence for specific reporter gene or GOI expression in Keratinocytes from granular layer	Inquiry
pPRO-KRT14 Plasmid	PTS048	Expression plasmid vector containing KRT14 promoter sequence for specific reporter gene or GOI expression in Keratinocytes from granular and basal layers	Inquiry
pPRO-ITGAM Plasmid	PTS049	Expression plasmid vector containing ITGAM promoter sequence for specific reporter gene or GOI expression in Monocytes, macrophages, granulocytes, natural killer cells	Inquiry
pPRO-UPK2 Plasmid	PTS050	Expression plasmid vector containing UPK2 promoter sequence for specific reporter gene or GOI expression in Bladder	Inquiry
pPRO-ENG Plasmid	PTS051	Expression plasmid vector containing ENG promoter sequence for specific reporter gene or GOI expression in Endothelial cells	Inquiry
pPRO-POU5F1 Plasmid	PTS052	Expression plasmid vector containing POU5F1 promoter sequence for specific reporter gene or GOI expression in Pluripotent cells (germ cells, ES cells, iPS cells)	Inquiry
pPRO-NANOG Plasmid	PTS053	Expression plasmid vector containing NANOG promoter sequence for specific reporter gene or GOI expression in Pluripotent cells (ES cells, iPS cells)	Inquiry

pPRO-SOX17 Plasmid	PTS054	Expression plasmid vector containing SOX17 promoter sequence for specific reporter gene or GOI expression in Endoderm	Inquiry
pPRO-ALB Plasmid	PTS055	Expression plasmid vector containing ALB promoter sequence for specific reporter gene or GOI expression in Liver	Inquiry
pPRO-CAMK2A Plasmid	PTS056	Expression plasmid vector containing CAMK2A promoter sequence for specific reporter gene or GOI expression in Brain	Inquiry
pPRO-TBG Plasmid	PTS057	Expression plasmid vector containing TBG promoter sequence for specific reporter gene or GOI expression in Liver	Inquiry
pPRO-TNNT2 Plasmid	PTS058	Expression plasmid vector containing TNNT2 promoter sequence for specific reporter gene or GOI expression in Heart	Inquiry
pPRO-CKM Plasmid	PTS059	Expression plasmid vector containing CKM promoter sequence for specific reporter gene or GOI expression in Heart	Inquiry
pPRO-ENO2 Plasmid	PTS060	Expression plasmid vector containing ENO2 promoter sequence for specific reporter gene or GOI expression in Brain	Inquiry
pPRO-GCG Plasmid	PTS061	Expression plasmid vector containing GCG promoter sequence for specific reporter gene or GOI expression in Brain and intestinal enteroendocrine cells but not islet cells	Inquiry